

MINI 3 Massage Device

Instructions

EN / CN / DE / FR / ES / IT / KR / PL / ID

- | | | |
|-------------------------------|----------------|------------------------------------|
| [1] Speed control button | [1] 调速键 | [1] Geschwindigkeitsregler |
| [2] Indicator | [2] 指示灯 | [2] Anzeigeleuchte |
| [3] Flat Head (Non-removable) | [3] 平头 (不可拆卸) | [3] Flacher Kopf (nicht abnehmbar) |
| [4] Cone head | [4] 圆锥头 | [4] Kegelförmiger Kopf |
| [5] U-shaped head | [5] U 型头 | [5] U-förmiger Kopf |
| [6] Head massaging head | [6] 头部按摩头 | [6] Kopf-Massagekopf |
| [7] Gentle head | [7] 轻柔按摩头 | [7] Sanfter Massagekopf |
| [8] USB-C cable | [8] Type-C 充电线 | [8] USB-C-Ladekabel |

- | | | |
|-------------------------------------|------------------------------------|--|
| [1] Bouton de réglage de la vitesse | [1] Botón de ajuste de velocidad | [1] Pulsante di regolazione della velocità |
| [2] Indicateur lumineux | [2] Indicador luminoso | [2] Indicatore luminoso |
| [3] Tête plate (non amovible) | [3] Cabezal plano (no desmontable) | [3] Testina piatta (non rimovibile) |
| [4] Tête conique | [4] Cabezal cónico | [4] Testina conica |
| [5] Tête en U | [5] Cabezal en forma de U | [5] Testina a forma di U |
| [6] Tête de massage pour la tête | [6] Cabezal para masaje de cabeza | [6] Testina massaggiante per la testa |
| [7] Tête de massage douce | [7] Cabezal para masaje suave | [7] Testina per massaggio delicato |
| [8] Câble de charge USB-C | [8] Cable de carga USB-C | [8] Cavo di ricarica USB-C |

- | | | |
|-------------------|-----------------------------------|--|
| [1] 속도 조절 버튼 | [1] Przycisk regulacji prędkości | [1] Tombol Pengatur Kecepatan |
| [2] 표시등 | [2] Wskaźnik LED | [2] Indikator Lampu |
| [3] 평면 헤드 (분리 불가) | [3] Głowica płaska (nierozłączna) | [3] Kepala Datar (Tidak Dapat Dilepas) |
| [4] 원뿔형 헤드 | [4] Głowica stożkowa | [4] Kepala Kerucut |
| [5] U 형 헤드 | [5] Głowica w kształcie U | [5] Kepala Bentuk U |
| [6] 두부 마사지 헤드 | [6] Głowica do masażu głowy | [6] Kepala Pijat untuk Kepala |
| [7] 부드러운 마사지 헤드 | [7] Głowica do delikatnego masażu | [7] Kepala Pijat Lembut |
| [8] Type-C 충전 케이블 | [8] Kabel ładowania USB-C | [8] Kabel Pengisian USB-C |

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Do not press the massage head firmly to block the motor rotation while using, otherwise it will reduce the service life of the product!

(1) Ten-minute protection: The machine will automatically enter standby mode after working continuously for 10 minutes. If the speed control button is slid within 3 minutes, the machine will resume operation and reset the timer. If there is no operation, the machine will automatically power off.

(2) Locked motor protection: If the massage head is obstructed and the rotation of the motor is blocked, the motor will stop, and the four gear indicator lights will flash three times to alert the user before automatically restarting. If blocking occurs five times within 5 minutes, the machine will automatically power off.

(3) Auto power off protection: The device will auto power off if there is no operation for 3 minutes.

※ In standby mode: The revolution speed is 0

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Specifications

Product name: KiCA Mini 3 Massage Device
Product model: KiCAM3
Battery capacity: 1100mAh
Operating temperature Range: 0°C ~+40°C
Charging time: About 2.5 Hours (5V --- 2A)
Charging time: About 7.5 Hours (Test in lab environment with Low speed)
Weight (Main body): 255±5g
Size (Main body): 89.2*35*138.1mm

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
X1

X1

This document is subject to change without notice.

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EN

Contraindications:

- | | |
|--|---|
| (1) Cancers | (7) Varicose veins |
| (2) Metal bone anastomosis (do not use the device directly in the fracture area) | (8) Acute inflammation |
| (3) Osteoporosis | (9) Do not massage the bony area |
| (4) Vein thrombosis and pulsation | (10) Do not massage surgical site |
| (5) Hypertension | (11) Do not massage the cervical vertebrae and spine |
| (6) Pregnant woman | (12) Kidney stones (do not use the device near the kidneys) |
| | (13) Do not massage one episode for more than 5 minutes |

Safety Instruction:

- (1) It is recommend to massage your body after wearing dry clothes or placing a clean and dry towel on massage part, then gently press and move the device.
- (2) High-speed and high-pressure massage can only be used to soft tissues without feeling pain or maladjustments. Do not massage head and any skeletal sites with high speed and high pressure.
- (3) Massage the same part of the body with high speed and high pressure for a long time will.
- (4) Please keep your fingers, hair or, any other body hair away from the extendable part of the device to prevent pinching and clipping.
- (5) Do not discard or misuse this device.
- (6) Do not modify this device by yourself without permission.
- (7) Do not operate or charge the device when nobody was around.
- (8) This device is for adults only.
- (9) People with contraindicated symptoms should use this device with doctor's permission.

Notice:

- (1) Make sure the massage head is not blocked by external force when the device is power on. If it is blocked, it will pause for a while and try to restart.
- (2) DO NOT contact the product with water or other liquid for the massage device is not waterproof.
- (3) DO NOT disassemble the product except marked detachable. It need send to after-sales or authorized service center to fix it if you accidentally disassemble it and cause abnormal work. The relevant costs are borne by users.
- (4) Prolonged continuous operation may cause the motor surface temperature to rise, please operate carefully.
- (5) DO NOT drop or strike the product. If the product is abnormal, contact After-sales support.

Storage and Maintenance:

- (1) Keep the product out of the reach of children and pets.
- (2) DO NOT leave the product near heat sources such as a furnace or heater. DO NOT leave the product inside of a vehicle on hot days.
- (3) Please storage the product in a dry environment.
- (4) Do not overcharge or overdischarge the battery, otherwise it will damage the battery cells. It is recommended to charge at least once a month when not in use.
- (5) Never use the product when the temperature is too high or too low.

Precautions when discarding the product:

The lithium ion battery used by the product is a recyclable and renewable valuable resource.

When discarding it, please have it handled by a professional agency that can recycle chargeable batteries. Before discarding the product, power it off and take the lithium ion battery out from it following the steps below. The battery is to be disposed of safely. Steps:

- (1) Disassemble the top housing with help of a professional tool.
- (2) Pull out the batter terminals from the PCBA.
- (3) Use a professional cutter to cut off the bottom housing and then take out the battery.

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